

Avalo Medical Trolley Series





University Hospitals Bristol NHS Foundation Trust

University Hospitals Bristol NHS Foundation Trust Resuscitation Services team applied for a major medical bid so they could replace and standardise all the resuscitation trolleys across their Trust.

Once granted a selection of trolleys were shortlisted and all hospital staff were invited to visit a large lecture theatre where the trolleys were displayed. This gave the staff from all clinical areas the opportunity to make a thorough assessment of the trolleys in terms of layout, drawer space, weight and manoeuvrability. After they had appraised all the options each person was able to vote for their preferred trolley and asked to give the reasons why. The Avalo Emergency trolley from Distinctive Medical was the most popular and won the vote by a large majority.

The key reasons why people chose the trolley were that it looks good, lightweight yet strong, easy to clean, great manoeuvrability yet stable, good drawer sizes and excellent flexibility in drawer configuration and layout. Staff also noted that they found the tamper-evident mechanism is easy to use.

Joanna Bruce-Jones the Resuscitation Lead at University Hospitals Bristol NHS Foundation Trust says the “trolleys have far exceeded our expectations”. While initially, the process of replacing the Trust’s resuscitation trolleys seemed like a daunting task with the assistance of Distinctive Medical it was made much simpler ensuring they got it right first time.

“After-sales service has been faultless with immediate responsiveness.

We could not recommend a better company to deal with or trolley to purchase.”

Joanna Bruce – Jones, Resuscitation Lead

How to order

If you would like to order or find out more about any Avalo medical trolleys our in-house team of specialists are on hand to help. Please call us on [01928 571801](tel:01928571801) or email info@distinctivemedical.com.

The Avalo Series Emergency Trolley Delivery Options

Standard Delivery

Your trolley will be delivered, securely packaged on a pallet, for the safe arrival to a loading bay location.

Packaging Removal*

All the benefits of the standard delivery option without the hassle of disposing of the packaging. We will take it back with us and dispose of it in an environmentally friendly way.

White Glove Delivery*

Our white glove delivery is second to none and by far the most efficient way to receive your trolley. We will fit any additional accessories onsite and deliver the completed trolley directly to the clinical area.

*Additional charges apply for this service



Call on **01928 571801**



Fax on **01928 571430**



Email on **info@distinctivemedical.com**



Live Chat on **www.distinctivemedical.com**

We accept payments by the following:



We also accept purchase orders from all organisations* as well as BACS payments.

* Subject to successful credit check

Avalo® Series

Resuscitation Emergency Trolley

The **Avalo Series Emergency Trolley** is one of the most acclaimed and trusted resuscitation trolleys on the UK market, with over 12,000 trolleys currently in use in hospitals up and down the country.

Equipped to support staff in an emergency with its light yet strong structure for the best manoeuvrability, through to the unique anti-tamper locking handle that allows for fast access to equipment making it one of the most innovative on the market.

The trolley is modular in design so can be configured and organised to match your hospital's requirements. Drawer dividers ensure each piece of equipment has its own place for easy location. Available in three sizes with a range of accessories to enhance user efficiency and support team performance.



Quick-Access lock offers the security you require and ensures quick access to emergency trolley supplies



Paediatric Emergency Trolley

The **Avalo Series Paediatric Trolley** has been designed and manufactured to the same high standards as the adult Avalo Resuscitation Trolley. Offering the same flexibility in design so it can be constructed to follow your hospital's procedure for treating paediatric patients in an emergency.

Available with the breakaway lock as standard giving immediate access to vital equipment during a rapid response - saving time when every moment matters.

The trolley comes decorated with child-friendly labels to help put patients at ease. All labels used are super durable and easy to decontaminate. This trolley is compatible with the accessories listed on page 13.



Difficult Airway Trolley

The **Avalo Difficult Airway Trolley** combines the manufacturing excellence of Avalo with the medical expertise of Distinctive Medical.

Designed exclusively for difficult airway situations, the trolley provides a standardised and structured way to store crucial equipment. The use of custom labels helps to guide clinicians through the process and identify the required equipment. Distinctive labelling means a Difficult Airway Trolley is instantly recognisable by both clinical and non-clinical staff.

The Difficult Airway Trolley comes with a bougie holder and difficult airway labels as standard.



Anaesthesia Trolley

The **Avalo Series Anaesthesia Trolley** is designed to meet the unique needs of your Anaesthetics Department. Boasting all the key features you would expect from an Avalo, but with the addition of an electronic keypad locking system, reducing the risk of a security breach.


Choose up to 5 different 4-digit codes which can be used for different members of staff or a simple key lock. Drug syringe label dispensers can also be attached to the accessory tower turning your Anaesthesia Trolley into a mobile, secure and efficient addition to your anaesthetics room.



Drawer Dividers



Removable Drawer Trays



Pullout Writing Surface



Isolation Trolley

The **Avalo Series Isolation Trolley** is made for the storage of personal protective equipment making it ideal for isolation wards. Available in a choice of drawer depths and optional drawer dividers for simple organisation of items.

Each drawer opens out all the way giving full access to all the contents of the trolley no matter where they are stored within it.

The outer surface of the trolley is made from a strong polymer plastic that can withstand continuous decontamination from strong chlorine-based cleaners.



Spread Caster Position provides stability control by moving the wheels out toward the corners of the cart frame

I.V. Therapy Trolley

The **Avalo Series I.V. Therapy Trolley** encourages the thorough organisation of items. Accessories options such as drawer dividers and storage bins allow for each element used for I.V. therapy to have its own place.

The top of the trolley provides a flat, spacious surface for the preparation of solutions. If extra space is needed an optional slide out panel can be added.



Integral Drawer Dividers offer the flexibility to organise and re-arrange drawer contents



Procedure Trolley

The **Avalo Series Procedure Trolley** is one of our most adaptable. With all the essential elements of a medical trolley it can be built up with combinations of accessories options to make it suitable for a specific department's needs.

Choose from three locking options. Breakaway lock – ideal if the trolley is intended for emergency use, keyless lock – keeps items protected and gives authorised members of staff instant access with a 4-digit code or the key lock – for advanced security.



Seamless Drawer Trays are durable and prevent issues associated with pop rivets or sharp edges found in metal drawers

Medical/Surgical Trolley

The **Avalo Series Med/Surg Trolley** has been designed to be compatible with the latest healthcare equipment and technology.

Available with a choice of add-ons so equipment can be stored or displayed in the most productive way to create an effective workflow. The efficacy of this trolley can be continued in the drawer space with a selection of drawer dividers and tray solutions to choose from.

Precision welded all-steel frame is virtually indestructible



Critical Care Trolley

The **Avalo Series Critical Care Trolley** has been designed in a way that appreciates the ever-changing demands within a hospital environment. Choose from a range of accessories options to help support staff with concise organisation.

Select different storage options to suit various sizes of equipment from small disposable items to large monitors. This trolley can be constructed to meet your exact specifications. Great for use in outreach environments.

See the full list of accessories on page 13 and see how the critical care trolley can work for you.

12 Year Limited Warranty on every Avalo Series Medical Trolleys



Avalo® Series Specifications

Integrated handle
on all trolleys

Custom drawer configurations
lets you easily configure and
re-configure drawers
to suit your requirements

12 year warranty
for peace of mind

Integrated bumpers
protects side panels
and castors



Available in 3 heights

- Tall: 110cm
- Standard: 100cm
- Compact: 91cm

Slide out surface

for even more working area


3 lock types available



Anti tamper seal lock



Key lock



Key pad

HDPE polymer

construction resistant to impacts and easily cleaned

Lockable caster

for exceptionally easy manoeuvring and securing



Built to last, designed to evolve

Our confidence in the reliability and durability of each Avalo trolley is reflected with a 12-year warranty. We understand that over the trolley's lifetime, your department's needs may change and for this reason, many elements are made to be re-configurable to evolve as you do.



Trolley Builder

Avalo has been a trusted brand in the UK market for over fifteen years.

Exceptionally high manufacturing quality, combined with intelligent design make the Avalo trolley a go-to choice for hospitals across the UK.

Available in a range of sizes, colours and with various accessories, Avalo trolleys can be modified to meet your exact requirements.

If we do not currently have an accessory that you require, please contact us.



Choose your height

Drawer Configuration*

Standard Height

Intermediate Height

Compact Height

Trolley Specifications

Applies to model types:

Emergency
Anaesthesia



	Surf. Height	Depth	Width
Standard Height:	110 cm	61 cm	79 cm
Intermediate Height:	100 cm	61 cm	79 cm
Compact Height:	91 cm	61 cm	79 cm

Applies to model types:

Procedure
Treatment
Critical Care
I.V. Therapy
Med/Surg



Drawer Specifications

Height	Depth	Width
7 cm	45 cm	55 cm
16 cm	45 cm	55 cm
25 cm	45 cm	55 cm

Applies to model type:

Isolation



*Standard drawer configurations shown, alternatives are available. Consult with your Distinctive Medical representative for bespoke cart handle configurations. Specifications, availability & components are subject to change without notice.

Choose your lock

All locking systems lock every drawer individually within 25mm of complete closure.



Anti-Tamper Seal Lock
Secure a trolley with a tamper evident seal for traceability. The seal is broken every time the trolley is opened by moving the handle down to unlock, allowing for quick access to equipment and supplies in emergency situations.



Keyless Access with Manual Relocking
An electronic keypad locking system, reducing the risk of a security breach. Choose up to 5 different 4-digit codes which can be used for different members of staff. Secure manually by pushing in the lock lever.



Key Access with Manual Relocking
The most simple but secure locking system by giving only those authorised access to the key. Ideal for trolleys that require advanced security. The removable lock system allows for flexibility by making changes to keys easy.

Choose your Colour



Red



Blue



Yellow



Green



Red / Grey



Blue / Grey



Yellow / Grey



Green / Grey



Grey / Light Grey

Choose your accessories



Glove box holders



Sharps/glove bin



Oxygen tank holder



Defibrillator shelf



I.V. Pole



Custom drawer labelling



Suction unit brackets



Writing surface




Electrical extension socket



Catheter holder



Tilting storage bins



Accessory bridge



Storage container



Label dispenser



Drawer dividers




Computer arm

Create your labels

We can create custom labels for trolleys that can help increase organisation, familiarity with contents and procedures or just simply brighten up its appearance.

Whatever your needs, our in-house Artwork and Printing team will be able to create the right labels for you.



Resuscitation Trolley



Difficult Airway Trolley



Paediatric Trolley

Talk to our team about your requirements on **01928 571801** or email info@distinctivemedical.com

Manage Your Compliance The Proactive Way

MyKitCheck is a proactive monitoring platform that puts compliance at your fingertips. Unlock all these benefits by taking your checks digital...



Digital Paperless Recording



Alerts For Expiring Equipment



Alerts For Missing Equipment



Digital Noticeboard



Automatic & Manual Online Ordering



Frequency of Checks Reduced



Real Time Compliance Monitoring



Advanced Tamper Proof Security



Instantly Locate Items for Recall



3 Seymour Court, Tudor Road, Manor Park, Runcorn, Cheshire WA7 1SY
01928 571 801 | info@distinctivemedical.com | www.distinctivemedical.com

Registered in the UK. Registered office: 3 Seymour Court, Tudor Road, Manor Park, Runcorn, Cheshire WA7 1SY
VAT Registration No. 775518887